




Product Datasheet

EIF5AL1 antibody (orb29960)

Description	Peptide affinity purified rabbit polyclonal antibody.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, IHC-P, WB
Immunogen	Synthetic Peptide
Target	EIF5AL1
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	16773
Uniprot ID	Q6IS14
NCBI	NP_001093162.1
Dilution Range	IHC-P: 1:10-50, WB: 1:1000, IF/ICC: 1:10-50
Expiration Date	12 months from date of receipt.



Western blot analysis of Ramos, A2058 ce...



Confocal immunofluorescence analysis of ...