

Product Datasheet

ALDH1L2 antibody (FITC) (orb29961)

Description	ALDH1L2 antibody (FITC)
Species/Host	Rabbit
Reactivity	Human
Conjugation	FITC
Immunogen	Recombinant Human Mitochondrial 10-formyltetrahydrofolate dehydrogenase protein (263-562AA)
Target	ALDH1L2
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q3SY69
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.