

Product Datasheet

TSGA10IP antibody (orb29976)

Description	TSGA10IP antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Recombinant Human Testis-specific protein 10-interacting protein (6-302AA)
Target	TSGA10IP
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q3SY00
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.