




## Product Datasheet

### PCDH1 antibody (orb29995)

|                            |                                                                                                                                                                              |
|----------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | Peptide affinity purified rabbit polyclonal antibody.                                                                                                                        |
| <b>Species/Host</b>        | Rabbit                                                                                                                                                                       |
| <b>Reactivity</b>          | Human                                                                                                                                                                        |
| <b>Conjugation</b>         | Unconjugated                                                                                                                                                                 |
| <b>Tested Applications</b> | IF, IHC-P, WB                                                                                                                                                                |
| <b>Immunogen</b>           | Synthetic Peptide                                                                                                                                                            |
| <b>Target</b>              | PCDH1                                                                                                                                                                        |
| <b>Form/Appearance</b>     | Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification. |
| <b>Note</b>                | For research use only                                                                                                                                                        |
| <b>Isotype</b>             | Rabbit IgG                                                                                                                                                                   |
| <b>Clonality</b>           | Polyclonal                                                                                                                                                                   |
| <b>MW</b>                  | 114743                                                                                                                                                                       |
| <b>Uniprot ID</b>          | <a href="#">Q08174</a>                                                                                                                                                       |
| <b>NCBI</b>                | <a href="#">NP_001265542.1</a> , <a href="#">NP_002578.2</a> , <a href="#">NP_115796.2</a> , <a href="#">NP_001265544.1</a>                                                  |
| <b>Dilution Range</b>      | IHC-P: 1:10-50, WB: 1:1000, IF/ICC: 1:10-50                                                                                                                                  |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                                                                                                              |



Western blot analysis of MDA-MB453, ZR-7...



Confocal immunofluorescence analysis of ...