

Product Datasheet

SLC25A52 antibody (orb30070)

Description	SLC25A52 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Recombinant Human Solute carrier family 25 member 52 protein (1-33AA)
Target	SLC25A52
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q3SY17
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.