



Product Datasheet

CYBA antibody (orb30124)

Description	Peptide affinity purified rabbit polyclonal antibody.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	Synthetic Peptide
Target	CYBA
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	21013
Uniprot ID	P13498
NCBI	NP_000092.2
Dilution Range	WB: 1:1000, IF/ICC: 1:10-50
Expiration Date	12 months from date of receipt.



Western blot analysis of CYBA antibody i...