




Product Datasheet

TCRB antibody (orb30172)

Description	Peptide affinity purified rabbit polyclonal antibody.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, WB
Immunogen	Synthetic Peptide
Target	TRBV7-9
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	13107
Uniprot ID	P04435
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.



WB analysis
of (lane 1)
Jurkat whole
cel...



WB analysis
of human
spleen
lysates usin...