

Product Datasheet

MOZ/MYST3 antibody (orb30176)

Description	Purified rabbit polyclonal antibody.
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	KAT6A
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is prepared by Saturated Ammonium Sulfate (SAS) precipitation followed by dialysis against PBS.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	225028
Uniprot ID	Q92794
NCBI	NP_006757.2 , NP_001092883.1 , NP_001092882.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

