

Product Datasheet

Integrin beta 1 antibody (orb301820)

Description	Mouse monoclonal antibody to Integrin beta 1
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IF, IHC-Fr
Immunogen	Please contact us for more information.
Preservatives	PBS with 0.02% sodium azide
Concentration	100 ug/ml
Storage	2-8°C
Note	For research use only
Isotype	IgG1,kappa
Clonality	Monoclonal
Dilution Range	„
Expiration Date	12 months from date of receipt.