

Product Datasheet

CD38 antibody (FITC) (orb301910)

Description	Mouse monoclonal antibody against CD38 conjugated to FITC
Species/Host	Mouse
Reactivity	Human
Conjugation	FITC
Tested Applications	FC
Immunogen	Please contact us for more information.
Preservatives	PBS with 4mg/ml rAlbumin and 0.1% sodium azide
Concentration	n/a
Storage	2-8°C
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.