

Product Datasheet

CD56 antibody (PE) (orb301998)

Description	Mouse monoclonal antibody against CD56 conjugated to PE
Species/Host	Mouse
Reactivity	Human
Conjugation	RPE
Tested Applications	FC
Immunogen	Please contact us for more information.
Preservatives	PBS with 4mg/ml rAlbumin, 10% sucrose and 0.1% sodium azide
Concentration	n/a
Storage	2-8°C
Note	For research use only
Isotype	IgG2a, k
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.