

Product Datasheet

Keratin 19 antibody (orb302050)

Description	Mouse monoclonal antibody to Keratin 19
Species/Host	Mouse
Reactivity	Human, Zebrafish
Conjugation	Unconjugated
Tested Applications	IHC-Fr, IHC-P, WB
Immunogen	Please contact us for more information.
Preservatives	PBS with 0.09% sodium azide
Concentration	n/a
Storage	2-8°C
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Dilution Range	,IHC-p (1:10-1:50),
Expiration Date	12 months from date of receipt.