

Product Datasheet

Epithelial Keratin type I antibody (orb302066)

Description	Mouse monoclonal antibody to Epithelial Keratin type I
Species/Host	Mouse
Reactivity	Bovine, Gallus, Human, Mouse, Rabbit, Rat
Conjugation	Unconjugated
Tested Applications	IHC-Fr, IHC-P, WB
Immunogen	Please contact us for more information.
Preservatives	PBS with 0.1% sodium azide
Concentration	1 mg/ml
Storage	2-8°C
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Dilution Range	IHC-P (1:150-1:200),IHC-Fr (1:150-1:200),WB (1:100)
Expiration Date	12 months from date of receipt.