

Product Datasheet

TLR3 antibody (orb302091)

Description	Mouse monoclonal antibody to TLR3
Species/Host	Mouse
Reactivity	Canine, Human, Mouse
Conjugation	Unconjugated
Tested Applications	FA, FC, IF, IHC-Fr, IHC-P, IP
Immunogen	Unknown or proprietary to Biorbyt and/or its suppliers
Preservatives	PBS with 0.1% rAlbumin and 0.02% sodium azide
Concentration	100 ug/ ml
Storage	2-8°C
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Dilution Range	
Expiration Date	12 months from date of receipt.