

Product Datasheet

Keratin 6 antibody (orb302093)

Description	Mouse monoclonal antibody to Keratin 6
Species/Host	Mouse
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IHC-Fr, IHC-P, WB
Immunogen	Please contact us for more information.
Preservatives	PBS with 0.5% rAlbumin and 0.09% sodium azide
Concentration	n/a
Storage	2-8°C
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Dilution Range	,IHC-P (1:10),
Expiration Date	12 months from date of receipt.