

Product Datasheet

Caldesmon antibody (orb302126)

Description	Mouse monoclonal antibody to Caldesmon
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	Please contact us for more information.
Preservatives	Tris Buffer, pH 7.3-7.7, containing 1% rAlbumin and 0.1% Sodium Azide
Concentration	n/a
Storage	2-8°C
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
Dilution Range	IHC-P (1:25-1:100)
Expiration Date	12 months from date of receipt.