



Product Datasheet

SLC36A3 antibody (orb30215)


Description	SLC36A3 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC
Immunogen	Recombinant Human Proton-coupled amino acid transporter 3 protein (1-46AA)
Target	SLC36A3
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q495N2
Dilution Range	IHC-P:1:20-200, IF:1:50-200
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of MCF-7 cel...



Immunohistochemical staining of human ki...



Immunohistochemical staining of human te...