

Product Datasheet

Cytomegalovirus antibody (orb302170)

Description	Mouse monoclonal antibody to Cytomegalovirus
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	Please contact us for more information.
Preservatives	Tris Buffer, pH 7.3-7.7, containing 1% rAlbumin and 0.1% Sodium Azide
Concentration	n/a
Storage	2-8°C
Note	For research use only
Isotype	IgG2a-k, IgG1-k
Clonality	Monoclonal
Dilution Range	IHC-P (1:10-1:20)
Expiration Date	12 months from date of receipt.