

Product Datasheet

ksp- Cadherin antibody (orb302193)

Description	Mouse monoclonal antibody to ksp-Cadherin
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	Please contact us for more information.
Preservatives	Tris Buffer, pH 7.3-7.7, containing 1% rAlbumin and 0.1% Sodium Azide
Concentration	n/a
Storage	2-8°C
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Dilution Range	IHC-P (1:100-1:500)
Expiration Date	12 months from date of receipt.