

Product Datasheet

Epstein Barr Virus antibody (orb302267)

Description	Mouse monoclonal antibody to Epstein Barr Virus
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	Please contact us for more information.
Preservatives	Tris Buffer, pH 7.3-7.7, containing 1% rAlbumin and 0.1% Sodium Azide
Concentration	n/a
Storage	2-8°C
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.