

Product Datasheet

Alpha-Tubulin antibody (orb302290)

Description	Mouse monoclonal antibody to Alpha-Tubulin
Species/Host	Mouse
Reactivity	Human, Mouse, Plant, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-Fr
Immunogen	Please contact us for more information.
Preservatives	PBS with 0.02% sodium azide
Concentration	100 ug/ml
Storage	2-8°C
Note	For research use only
Isotype	IgM
Clonality	Monoclonal
Dilution Range	,
Expiration Date	12 months from date of receipt.