

Product Datasheet

TLR2 antibody (Biotin) (orb302298)

Description	Mouse monoclonal antibody against TLR2 conjugated to Biotin
Species/Host	Mouse
Conjugation	Biotin
Immunogen	Unknown or proprietary to Biorbyt and/or its suppliers
Preservatives	PBS with 0.1% rAlbumin and 0.02% sodium azide
Concentration	100 ug/ ml
Storage	2-8°C
Note	For research use only
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.