



Product Datasheet

PCDH11X antibody (orb3023)

Description	PCDH11X antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	KLH conjugated synthetic peptide derived from human PCDH11X (601-700/1347 aa)
Target	PCDH11X
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	145 kDa
Uniprot ID	Q9BZA8, Q9BZA7
Dilution Range	WB=1:500-2000
Expiration Date	12 months from date of receipt.



Western blot: Lane
1: Mouse
Hippocampus; La...



Western blot: Lane
1: Human SH-SY5Y
cell; L...