

Product Datasheet

AFP antibody (orb302312)

Description	Mouse monoclonal antibody to AFP
Species/Host	Mouse
Reactivity	Canine, Human, Monkey, Porcine
Conjugation	Unconjugated
Tested Applications	ELISA, FC, IF, IHC-Fr, IHC-P
Immunogen	Please contact us for more information.
Preservatives	PBS with 0.02% sodium azide
Concentration	100 ug/ml
Storage	2-8°C
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
Dilution Range	
Expiration Date	12 months from date of receipt.