

Product Datasheet

CD62I antibody (PE) (orb302320)

Description	Mouse monoclonal antibody against CD62I conjugated to PE
Species/Host	Mouse
Reactivity	Bovine, Canine, Human, Primate
Conjugation	RPE
Tested Applications	FC
Immunogen	Please contact us for more information.
Preservatives	PBS with 0.09% sodium azide and 1% rAlbumin, 5% sucrose
Concentration	n/a
Storage	2-8°C
Note	For research use only
Isotype	IgG2b
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.