

Product Datasheet

PEnicillin antibody (orb302321)

Description	Mouse monoclonal antibody to PEnicillin
Species/Host	Mouse
Reactivity	Other
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Please contact us for more information.
Preservatives	PBS with 0.02% sodium azide
Concentration	100 ug/ml
Storage	2-8°C
Note	For research use only
Isotype	IgG1-kappa
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.