

## Product Datasheet

### CD66aCD antibody (orb302330)

|                            |                                         |
|----------------------------|-----------------------------------------|
| <b>Description</b>         | Mouse monoclonal antibody to CD66aCD    |
| <b>Species/Host</b>        | Rat                                     |
| <b>Reactivity</b>          | Human, Primate                          |
| <b>Conjugation</b>         | Unconjugated                            |
| <b>Tested Applications</b> | FC, IHC-Fr, IHC-P                       |
| <b>Immunogen</b>           | Please contact us for more information. |
| <b>Preservatives</b>       | PBS with 0.09% sodium azide             |
| <b>Concentration</b>       | 1.0 mg/ml                               |
| <b>Storage</b>             | 2-8°C                                   |
| <b>Note</b>                | For research use only                   |
| <b>Isotype</b>             | IgG2a                                   |
| <b>Clonality</b>           | Monoclonal                              |
| <b>Dilution Range</b>      | „                                       |
| <b>Expiration Date</b>     | 12 months from date of receipt.         |