

Product Datasheet

CD103 antibody (PE) (orb302341)

Description	Mouse monoclonal antibody against CD103 conjugated to PE
Species/Host	Mouse
Reactivity	Human
Conjugation	RPE
Tested Applications	FC
Immunogen	Please contact us for more information.
Preservatives	PBS with 0.09% sodium azide and 1% rAlbumin, 5% sucrose
Concentration	n/a
Storage	2-8°C
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.