

Product Datasheet

Gremlin 1 antibody (orb30236)

Description	Rabbit polyclonal antibody to Gremlin 1
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	GREM1
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	20697
Uniprot ID	O60565
NCBI	NP_037504.1 , NP_001178252.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

