

Product Datasheet

CD88 antibody (orb302379)

Description	Mouse monoclonal antibody to CD88
Species/Host	Mouse
Reactivity	Human, Primate
Conjugation	Unconjugated
Tested Applications	FC, IHC-Fr, IHC-P, WB
Immunogen	Please contact us for more information.
Preservatives	PBS with 0.09% sodium azide
Concentration	1 mg/ml
Storage	2-8°C
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
Dilution Range	...
Expiration Date	12 months from date of receipt.