

Product Datasheet

CD226 antibody (orb302408)

Description	Mouse monoclonal antibody to CD226
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IP
Immunogen	Please contact us for more information.
Preservatives	PBS with 0.09% sodium azide
Concentration	1.0 mg/ml
Storage	2-8°C
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Dilution Range	,
Expiration Date	12 months from date of receipt.