

Product Datasheet

Mouse IgG1 antibody (Biotin) (orb302474)

Description	Mouse monoclonal antibody against Mouse IgG1 conjugated to Biotin
Species/Host	Rat
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-Fr
Immunogen	Please contact us for more information.
Preservatives	PBS with 0.1% sodium azide
Concentration	1 mg/ml
Storage	-20°C
Note	For research use only
Isotype	IgGk
Clonality	Monoclonal
Dilution Range	,
Expiration Date	12 months from date of receipt.