

## Product Datasheet

### Mouse IgG2a antibody (Peroxidase) (orb302499)

|                            |                                                                           |
|----------------------------|---------------------------------------------------------------------------|
| <b>Description</b>         | Mouse monoclonal antibody against<br>Mouse IgG2a conjugated to Peroxidase |
| <b>Species/Host</b>        | Rat                                                                       |
| <b>Reactivity</b>          | Mouse                                                                     |
| <b>Conjugation</b>         | Unconjugated                                                              |
| <b>Tested Applications</b> | ELISA, IHC-Fr                                                             |
| <b>Immunogen</b>           | Please contact us for more information.                                   |
| <b>Preservatives</b>       | PBS pH 7.2 / 50% (v/v) Glycerol                                           |
| <b>Concentration</b>       | 1 mg/ml                                                                   |
| <b>Storage</b>             | -20°C                                                                     |
| <b>Note</b>                | For research use only                                                     |
| <b>Isotype</b>             | IgG1k                                                                     |
| <b>Clonality</b>           | Monoclonal                                                                |
| <b>Dilution Range</b>      | ,                                                                         |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                           |