

Product Datasheet

Dinitrophenyl Hapten antibody (orb302500)

Description	Mouse monoclonal antibody to Dinitrophenyl Hapten
Species/Host	Rat
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Please contact us for more information.
Preservatives	PBS with 0.1% sodium azide
Concentration	1 mg/ml
Storage	-20°C
Note	For research use only
Isotype	IgG1k
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.