

Product Datasheet

IL1 Receptor antibody (orb302506)

Description	Mouse monoclonal antibody against IL1 Receptor conjugated to phospho-80
Species/Host	Rat
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IP, WB
Immunogen	Please contact us for more information.
Preservatives	PBS with 0.1% rAlbumin and 0.02% sodium azide
Concentration	100 ug/ ml
Storage	2-8°C
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
Dilution Range	„
Expiration Date	12 months from date of receipt.