

Product Datasheet

IFN-Alpha II antibody (orb302519)

Description	Mouse monoclonal antibody to IFN-Alpha II
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-Fr, WB
Immunogen	Please contact us for more information.
Preservatives	PBS with 0.02% sodium azide
Concentration	100 ug/ml
Storage	2-8°C
Note	For research use only
Isotype	IgG1,kappa
Clonality	Monoclonal
Dilution Range	„
Expiration Date	12 months from date of receipt.