

Product Datasheet

CD90 antibody (orb302529)

Description	Mouse monoclonal antibody to CD90
Species/Host	Mouse
Reactivity	Human, Primate
Conjugation	RPE
Tested Applications	FC
Immunogen	Please contact us for more information.
Preservatives	PBS with 0.09% sodium azide and 1% rAlbumin, 5%sucrose
Concentration	n/a
Storage	2-8°C
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.