

Product Datasheet

Liver Fatty Acid Binding Protein antibody (orb302563)

Description	Mouse monoclonal antibody to Liver Fatty Acid Binding Protein
Species/Host	Mouse
Reactivity	Canine, Human, Porcine, Primate, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-Fr, IHC-P, WB
Immunogen	Please contact us for more information.
Preservatives	PBS with 0.1% rAlbumin and 0.02% sodium azide
Concentration	100 ug/ ml
Storage	2-8°C
Note	For research use only
Isotype	IgG2b
Clonality	Monoclonal
Dilution Range	...
Expiration Date	12 months from date of receipt.