

Product Datasheet

MUC1 (polymorphic epithelial mucin) antibody (orb302584)

Description	Mouse monoclonal antibody to MUC1 (polymorphic epithelial mucin)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, FC, IF, IHC-Fr, IHC-P
Immunogen	Please contact us for more information.
Preservatives	PBS with 0.02% sodium azide
Concentration	100 ug/ml
Storage	2-8°C
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Dilution Range	,,,,
Expiration Date	12 months from date of receipt.