

Product Datasheet

MUC5AC antibody (orb302595)

Description	Mouse monoclonal antibody to MUC5AC
Species/Host	Mouse
Reactivity	Bovine, Feline, Human, Monkey, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-Fr, IHC-P, WB
Immunogen	Please contact us for more information.
Preservatives	PBS with 0.02% sodium azide
Concentration	100 ug/ml
Storage	2-8°C
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Dilution Range	...
Expiration Date	12 months from date of receipt.