

Product Datasheet

MMP-3 antibody (orb302597)

Description	Mouse monoclonal antibody to MMP-3
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-Fr, IHC-P, WB
Immunogen	Please contact us for more information.
Preservatives	Culture supernatant with 0.05% sodium azide
Concentration	n/a
Storage	2-8°C
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Dilution Range	„
Expiration Date	12 months from date of receipt.