

Product Datasheet

T-cell Receptor antibody (orb302603)

Description	Mouse monoclonal antibody to T-cell Receptor
Species/Host	Mouse
Reactivity	Rat
Conjugation	Unconjugated
Tested Applications	FC, IHC-Fr
Immunogen	Please contact us for more information.
Preservatives	Culture supernatant with 0.05% sodium azide
Concentration	n/a
Storage	2-8°C
Note	For research use only
Isotype	IgG2b
Clonality	Monoclonal
Dilution Range	,
Expiration Date	12 months from date of receipt.