

Product Datasheet

R73 & CC52 antibody (orb302606)

Description	Mouse monoclonal antibody to R73 & CC52
Species/Host	Mouse
Reactivity	Rat
Conjugation	Unconjugated
Tested Applications	FA
Immunogen	Please contact us for more information.
Preservatives	Culture supernatant with 0.05% sodium azide
Concentration	n/a
Storage	2-8°C
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.