

Product Datasheet

Ep-Cam antibody (orb302636)

Description	Mouse monoclonal antibody to Ep-Cam
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, FC, IHC-P, WB
Immunogen	Please contact us for more information.
Preservatives	PBS with 15 mM sodium azide, pH 7.4
Concentration	250 ug/ ml
Storage	2-8°C
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Dilution Range	...
Expiration Date	12 months from date of receipt.