

Product Datasheet

SIGN-R1 antibody (orb302662)

Description	Mouse monoclonal antibody to SIGN-R1
Species/Host	Rat
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	FA, FC, IHC-Fr
Immunogen	Please contact us for more information.
Preservatives	Sterile cell culture medium with 0.02% sodium azide
Concentration	200 ug/ ml
Storage	2-8°C
Note	For research use only
Isotype	IgM
Clonality	Monoclonal
Dilution Range	„
Expiration Date	12 months from date of receipt.