

Product Datasheet

LTA antibody (orb302669)

Description	Mouse monoclonal antibody to LTA
Species/Host	Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, FC, IF, IHC-Fr, WB
Immunogen	Please contact us for more information.
Preservatives	PBS with 0.1% rAlbumin and 0.02% sodium azide
Concentration	> 200 ug/ ml
Storage	2-8°C
Note	For research use only
Isotype	IgG3
Clonality	Monoclonal
Dilution Range	,,,,
Expiration Date	12 months from date of receipt.