

Product Datasheet

HSV1 + HSV2 antibody (orb302686)

Description	Mouse monoclonal antibody to HSV1 + HSV2
Species/Host	Mouse
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA, ICC
Immunogen	Please contact us for more information.
Preservatives	Phosphate buffered saline (PBS) solution with 15 mM sodium azide
Concentration	1 mg/ml
Storage	2-8°C
Note	For research use only
Isotype	IgG2b
Clonality	Monoclonal
Dilution Range	,
Expiration Date	12 months from date of receipt.