

Product Datasheet

Cellular Fibronectin antibody (orb302696)

Description	Mouse monoclonal antibody to Cellular Fibronectin
Species/Host	Mouse
Reactivity	Gallus, Guinea pig, Human, Rabbit, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-Fr, WB
Immunogen	Please contact us for more information.
Preservatives	PBS with 1% rAlbumin & 0.1% sodium azide
Concentration	100 ug/ ml
Storage	2-8°C
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Dilution Range	„
Expiration Date	12 months from date of receipt.