

Product Datasheet

Bromo-deoxy-Uridine antibody (orb302707)

Description	Mouse monoclonal antibody to Bromo-deoxy-Uridine
Species/Host	Mouse
Reactivity	Other
Conjugation	Unconjugated
Tested Applications	FC, ICC, IHC, IHC-Fr, IHC-P
Immunogen	Please contact us for more information.
Preservatives	PBS with 0.09% sodium azide
Concentration	1 mg/ml
Storage	2-8°C
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Dilution Range	FC (1:100-1:200),,,,
Expiration Date	12 months from date of receipt.