

Product Datasheet

Ribonucleoprotein Particles antibody (orb302724)

Description	Mouse monoclonal antibody to Ribonucleoprotein Particles
Species/Host	Mouse
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	FC, IF, IHC-Fr, IHC-P
Immunogen	Please contact us for more information.
Preservatives	PBS with 0.02% sodium azide
Concentration	100 ug/ml
Storage	2-8°C
Note	For research use only
Isotype	IgM
Clonality	Monoclonal
Dilution Range	...
Expiration Date	12 months from date of receipt.